
# TestPack+Plus with OBC Strep A

**[www.alere.com](http://www.alere.com)**

## INTENDED USE

**Alere™ TestPack Plus Strep A** with On Board Controls (OBC)

(**Alere™ TestPack Strep A**) is a rapid immunoassay for the qualitative detection of Group A Streptococcal (Group A Strep) antigen from throat swab specimens as an aid to the diagnosis of Group A Strep pharyngitis or for culture confirmation. For professional *in vitro* diagnostic use only.

## INTRODUCTION

Beta-haemolytic Group A *Streptococcus* is a major cause of upper respiratory infections in humans. The most commonly occurring Group A Streptococcal disease is pharyngitis. The symptoms of this, if left untreated, can become more severe and further complications such as acute rheumatic fever, toxic shock-like syndrome and glomerulonephritis can develop<sup>1</sup>. Rapid identification can facilitate clinical management to prevent disease progression.

**Alere™ TestPack Strep A** employs the Lancefield grouping method of identification, where Streptococcal groups are identified according to their cell wall antigens, which are species specific<sup>2,3</sup>.

Conventional methods used to identify Group A *Streptococcus* involve the isolation and subsequent identification of the organisms, which can take 24-48 hours to complete<sup>2,4</sup>.

**Alere™ TestPack Strep A** detects Group A Streptococci directly from throat swabs so that more rapid results are achieved. The test detects bacterial antigen from swabs, therefore it is possible to detect Group A *Streptococcus*, which may fail to grow in culture.

## TEST PRINCIPLE

The Streptococcal Group A specific antigen is extracted from the throat swab using Reagent 1 and Reagent 2. Following this, Reagent 3 is added to neutralise the acid formed by Reagents 1 and 2. The extract is then dropped into the Sample Well of the reaction disc where it mobilises the Group A Strep

antibody-coated colloid, which migrates through the membrane to the End of Assay Window.

If Group A Streptococcal antigen is present in the extract it will form a complex with the antibody-colloid. The antibody-colloid complex migrates through the membrane and is then captured by the Group A Strep antibody in the result window, providing a visual indication of the presence of antigen.

A pink/red Plus Sign "+" appearing in the Result Window indicates the presence of Group A Strep antigen. A Minus Sign "-" indicates no antigen was detected.

## KIT CONTENTS AND STORAGE


**REAGENT 1** 2.0 M Sodium Nitrite (& Xylenol Orange) (5 mL/10 mL)

**REAGENT 2** 1.0 M Acetic Acid (5 mL/10 mL)

**REAGENT 3** 1.0 M Tris Buffer (Preservative: Sodium Azide) (5 mL/10 mL)

- 20/40 reaction discs
- 20/40 extraction tubes
- 20/40 sterile Polyester (Dacron) tipped swabs
- One package insert

Store at 2-30°C. Do not use after the expiry date.

## PRECAUTIONS

Standard guidelines for handling infectious agents should be observed throughout all procedures.

1. It is recommended that disposable gloves be worn whilst handling specimens and reagents.
2. Do not mix extraction reagents and reaction discs from different lots.
3. Do not mix reagent bottle caps.
4. Reagent 2 or Reagent 1 combined with Reagent 2 is acidic. Avoid contact with eyes or mucous membranes. In the event of accidental contact, wash thoroughly with water.
5. Reagent 3 contains sodium azide, which on contact with lead and copper plumbing may react to form explosive metal azides. Use large volumes of water to flush reagents on disposal.
6. Properly dispose of all contaminated waste such as swabs, reaction discs and extraction tube.
7. Do not use reaction discs or swabs that have become wet or if the pouch has been opened or damaged.
8. Recap the reagent bottles between uses.

Hazard information for the components under applicable European Community (EC) Directives is as follows:

**REAGENT 1** R25, S1/2-45-60

Toxic: Contains sodium nitrite.

R25 Toxic if swallowed.

S1/2-45-60 Keep locked up and out of the reach of children. In case of accident, or if you feel unwell, seek medical advice immediately (show the label where possible). This material and its container must be disposed of as hazardous waste.

**REAGENT 3** R22, S2-35-46

Harmful: Contains sodium azide

R22 Harmful if swallowed.

S2-35-46 Keep out of the reach of children. This material and its container must be disposed of in a safe way. If swallowed, seek medical advice immediately and show this container or label.

Safety data sheet available for professional user on request.

#### SPECIMEN COLLECTION & STORAGE

A specimen should be obtained by standard throat swab collection Methods<sup>5,6</sup>. Use Polyester (Dacron) tipped swabs only. Do not use calcium alginate swabs, semi-solid transport media, or media containing charcoal.

If swabs are not to be processed immediately, they can be stored in clean, dry sealable plastic tubes refrigerated (2-8°C) for up to 72 hours before testing.

#### PROCEDURE

If the reagents, reaction discs and patient specimen are stored refrigerated allow them to equilibrate to room temperature (18-30°C) for a minimum of 30 minutes before beginning the assay. Do not open foil pouches until ready to perform the assay.

#### Extraction

All drops must be free falling with the reagent bottles held vertically.

1. Add 3 drops of Reagent 1 to an extraction tube. This solution must be pink.
2. Add 3 drops of Reagent 2 to the same tube. This solution must turn yellow.
3. Place the specimen swab into the tube. Twirl swab and mix well. Let stand for a minimum of 1 minute, but no longer than 30 minutes.
4. Add 3 drops of Reagent 3 to the same tube. Twirl swab and mix well. This solution must turn pink. Occasionally it may be necessary to add a fourth drop of Reagent 3 to obtain colour change.
5. **Thoroughly express the liquid from the swab by squeezing the tube with the thumb and forefinger, then rotating it as it is withdrawn.**
6. Discard the swab safely.
7. Apply the dropper tip to the extraction tube.

The extract is stable for up to 72 hours refrigerated (2-8°C) in a covered tube.

#### Test Procedure

Remove the reaction disc from its foil pouch. Label with patient or control identifications. Place on a clean, flat, dry surface.

8. Add the entire contents of the extraction tube drop-wise into the Sample Well on the reaction disc.
9. Wait for a pink/red colour in the End of Assay Window to appear (approximately 5 minutes). Read results. **Do not interpret results after 10 minutes from addition of extract.**

NOTE: To use Alere™ TestPack Strep A as a culture confirmation test, remove an isolated beta-haemolytic colony from the culture plate with a clean Polyester (Dacron) tipped swab and subject it to the complete test procedure, above.

#### INTERPRETATION OF RESULTS


##### On Board Controls

Alere™ TestPack Strep A utilises five On Board Controls to ensure the assay is functioning properly.

- The Extraction Reagent Control is demonstrated by colour changes during the swab extraction process to indicate that Reagents 1, 2 and 3 are added in the correct order. Reagent 1 is pink, but will change to yellow when Reagent 2 is added. After adding Reagent 3 the solution changes back from yellow to pink. If either of the colour changes does not occur, the test is invalid.


- The Positive On Board Control (POS CTL "✓") indicates that both the antibody-colloid complex and capture antibody systems are functional. The POS CTL "✓" must appear for the result to be valid.

- The Negative On Board Control (NEG CTL "X") indicates that the test specimen may contain a non-specific entity that could cause a false positive result. If the NEG CTL "X" appears in the Result Window, the result is invalid.

- The Minus Sign “–” indicates that migration of the extract has occurred. The Minus Sign “–” must appear for the result to be valid.

- The End of Assay Window must change to a red/pink colour, which indicates that the test is complete and ready to read. The End of Assay colour change must occur for the result to be valid. **Do not interpret results after 10 minutes from addition of extract.**

### Positive Result

A positive result consists of one vertical line (Patient Bar) and one horizontal line (Minus Sign "-") in the Result Window forming a Plus Sign "+". Pink/red colour (darker than the background) on the Patient Bar is interpreted as a positive result even if it has less colour than the Minus Sign. Randomly occurring red dots should not be evaluated in the interpretation of results.

### Negative Result

A negative result is indicated by a horizontal line (Minus Sign "-") in the Result Window. A negative result means that no Group A Strep antigen was detected, or that the levels of antigen in the specimen were below the detection limit of the assay.

If the test is invalid, retest with a new specimen and a new reaction disc, ensuring sufficient specimen addition and read for presence of the above controls. If the problem persists, call your local distributor.

### EXTERNAL QUALITY CONTROL

Good laboratory practice recommends the use of control materials to ensure proper kit performance. Each laboratory should refer to guidelines established internally and by local, national or other accrediting organisations.

Alternatively, use a sterile swab to obtain a confirmed Group A Streptococcal colony from a culture plate and subject this to the complete extraction and test procedure. A fresh sterile swab can be used as a negative control. If the external controls do not produce the expected result, the test is invalid and the patient test result should not be reported. Repeat testing of Positive/Negative Controls and patient specimens with new reaction discs. If the problem persists, call your local distributor.

### LIMITATIONS OF THE TEST

1. The use of swab types other than Polyester (Dacron) tipped swabs, swabs taken from sites other than the back of the throat or the use of other specimens such as saliva, sputum, or urine has not been established.

2. This test does not differentiate between carrier and infected individuals. Pharyngitis may be caused by organisms other than Group A *Streptococcus*. Performance characteristics in populations other than those studied during the clinical investigation have not been determined.
3. A negative result may be obtained if the amount of extracted antigen is below the sensitivity of the test.
4. False negative results may occur from incorrectly taken/extracted specimens.
5. Additional follow-up testing using the culture method is required if the result is negative and clinical symptoms persist.

### EXPECTED VALUES

It is believed that approximately 19% of all upper respiratory tract infections are caused by Group A *Streptococcus*<sup>7</sup>. Group A Strep associated pharyngitis displays a seasonal variation and is most prevalent during the winter and early spring. Certain populations are at higher risk of infection for example in schools, nursing homes and hospitals, and clustering of cases does occur<sup>8,9</sup>.

### CALIBRATION

**Alere™ TestPack Strep A** is calibrated using in-house standards produced from dilutions of a Group A Streptococcal antigen.

### PERFORMANCE CHARACTERISTICS

#### Clinical performance of Alere™ TestPack Strep A compared with standard sheep blood agar culture

In a multi-centre field evaluation, throat swabs were collected from 369 children and adults presenting with symptoms of pharyngitis. All swabs were tested on the day of collection and were used to inoculate a sheep blood agar (SBA) plate, prior to being tested using **Alere™ TestPack Strep A**. Plates were incubated for 24-48 hours at 35°C with 5-10% CO<sub>2</sub>. Presumptive group A *Streptococcus* colonies on SBA culture plates were confirmed using a commercially available Streptococcal latex grouping test.

The results are summarised below:

### Individual test results compared to standard SBA culture:

	SBA +	SBA -	Total
<b>TestPack Strep A +</b>	122	4	126
<b>TestPack Strep A -</b>	3	240	243
<b>Total</b>	125	244	369

Sensitivity: 122/125 = 97.6% (95% confidence interval [CI]: 93.1-99.5%)

Specificity: 240/244 = 98.4% (95% confidence interval [CI]: 95.9-99.6%).

### Individual test results compared to SBA culture density:

For visually positive SBA plates, clinicians recorded the density and growth of presumptive group A *Streptococcus* colonies. The results from SBA culture and the corresponding **Alere™ TestPack Strep A** rapid results are compared below:

Culture Density	TestPack Strep A +
1+	82% (9/11)
2+	96% (23/24)
3+	100% (39/39)
4+	100% (51/51)

## CROSS REACTIVITY

No cross reactivity was found when Alere™ TestPack Strep A was tested with the bacteria listed below. Organisms were tested at  $1 \times 10^8$  organisms per ml, with the exception of *Staphylococcus aureus*, which was tested at  $1 \times 10^9$  organisms per ml.

*Arcanobacterium haemolyticum*  
*Bordetella pertussis*  
*Candida albicans*  
*Corynebacterium diphtheriae*  
*Escherichia coli*  
*Fusobacterium necrophorum*  
*Haemophilus influenzae*  
*Haemophilus parahaemolyticus*  
*Klebsiella pneumoniae*  
*Moraxella catarrhalis*  
*Moraxella lacunata*  
*Neisseria gonorrhoeae*  
*Neisseria lactamica*  
*Neisseria meningitidis*  
*Neisseria sicca*  
*Neisseria subflava*  
*Proteus vulgaris*  
*Pseudomonas aeruginosa*  
*Serratia marcescens*  
*Staphylococcus aureus*  
*Staphylococcus aureus* (Cowan's serotype 1)  
*Staphylococcus epidermidis*  
*Staphylococcus haemolyticus*  
*Staphylococcus saprophyticus*  
*Streptococcus Groups B, C, D, F, G*  
*Streptococcus mitis*  
*Streptococcus mutans*  
*Streptococcus oralis*  
*Streptococcus pneumoniae*  
*Streptococcus salivarius*

*Streptococcus sanguis*  
*Yersina enterocolitica*

## ADVICE LINE

For further information, please contact your distributor, or call Alere Technical Specialists:

**Africa, Russia & CIS:** Tel: +972 8 9429 683  
Email: ARCSproductsupport@alere.com

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The Alere Logo and Alere are trademarks of the Alere group of companies.
© 2011 Alere. All rights reserved.
Sterile Polyester (Dacron) tipped swabs manufactured for Alere Medical Co., Ltd.
Polyesterummantelte (Dacron) Abstrichtupfer, hergestellt von Alere Medical Co., Ltd.
Sterile pødepinde med polyesterspids (Dacron) fremstillet for Alere Medical Co., Ltd.
Hisopos estériles con punta de poliéster (Dacron) fabricados por Alere Medical Co., Ltd.