


TestPack[®] Plus with OBC hCG Urine

www.alere.com

INTENDED USE

Alere™ TestPack Plus hCG Urine with On Board Controls (OBC) (Alere™ TestPack hCG Urine) is a rapid immunoassay for the qualitative detection of human chorionic gonadotropin (hCG) in urine for early detection of pregnancy. For professional *in vitro* diagnostic use only.

INTRODUCTION

Human chorionic gonadotropin (hCG) is a glycoprotein hormone that is produced by the blastocyst.¹

The background concentration of hCG in urine increases with age, but is normally <5mIU/ml in women of childbearing age². This rapidly increases after conception, reaching 50-250mIU/ml by the day of the expected period and peaks at approximately 100,000 to 200,000mIU/ml during the first trimester^{3,4}. The sudden rapid rise in concentration of hCG in urine following conception makes it an excellent marker for pregnancy.

TEST PRINCIPLE

Urine is added to the Sample Well of the reaction disc with the aid of a transfer pipette, and allowed to migrate through the membrane until it reaches the End of Assay Window. As the urine proceeds through the membrane, it mobilises the anti-hCG monoclonal antibody-colloid. If hCG is present in the specimen, it will form a complex with the antibody-colloid. The antibody-colloid complex migrates through the membrane and is captured by the immobilised anti-hCG monoclonal antibody in the result window, providing a visual indication of the presence of hCG.


The result must be read when the End of Assay Window has turned pink/red. If hCG is present in the urine at levels of 25mIU/ml or greater, a Plus Sign "+" appears in the Result Window. A Minus Sign "-" indicates no hCG was detected.

KIT CONTENTS AND STORAGE

- 20 reaction discs
- Pack of 20 transfer pipettes
- One package insert

Store at 2-30°C. Do not use after the expiry date.

Safety data sheet available for professional user on request.

**PRECAUTIONS**

Standard guidelines for handling infectious agents should be observed throughout all procedures.

1. Do not use reaction discs that have become wet or if the foil pouch has been opened or damaged.
2. Properly dispose of all contaminated waste such as reaction discs and transfer pipettes.

SPECIMEN COLLECTION & STORAGE

A urine specimen collected at any time of the day is suitable, but a first morning specimen is recommended for early detection of pregnancy, as it should contain the highest concentration of hCG⁵.

Urine specimens must be collected in **clean**, dry, plastic or glass containers. Specimens may be stored in the refrigerator (2-8°C) for up to 48 hours, or frozen **once** (-20°C) for up to 3 months. Specimens should not be repeatedly frozen and thawed.

Alere™ TestPack hCG Urine has not been validated for use with specimens containing preservatives other than sodium azide (0.1%).

No centrifugation or filtration of specimens is required prior to testing. However, particulate matter in specimens should be allowed to settle before removing a sediment-free aliquot for testing.

TEST PROCEDURE

Allow the reaction disc and patient specimen to equilibrate to room temperature (18-30°C) for a minimum of 30 minutes before beginning the assay. Do not open foil pouches until ready to perform the assay.

1. Remove the reaction disc from its foil pouch. Label with patient or control identifications. Place on a clean, flat, dry surface.
2. Draw specimen to the line marked on the transfer pipette. Dispense the entire contents drop wise into the Sample Well on the reaction disc.
3. Wait for a pink/red colour in the End of Assay Window to appear (approximately 5 minutes). Read results. **Ignore any results after this time.**


INTERPRETATION OF RESULTS

On Board Controls

Alere™ TestPack hCG Urine utilises four On Board Controls to ensure the assay is functioning properly.

- The **Positive On Board Control** (POS CTL "✓") indicates that both the antibody-colloid complex and capture antibody systems are functional. The POS CTL "✓" must appear for the result to be valid.

- The **Negative On Board Control** (NEG CTL "X") indicates that the test specimen may contain a non-specific entity that could cause a false positive result. If the NEG CTL "X" appears in the Result Window, the result is invalid.

- The **Minus Sign** “-” indicates that migration of the specimen has occurred. The Minus Sign “-” must appear for the result to be valid.
- The **End of Assay Window** must change to a red/pink colour, which indicates that the test is complete and ready to read. The End of Assay colour change must occur for the result to be valid.

Positive Result

A positive result consists of one vertical line (Patient Bar) and one horizontal line (Minus Sign) in the Result Window forming a Plus Sign "+". Pink/red colour (darker than the background) on the Patient Bar is interpreted as a positive result even if it has less colour than the Minus Sign. Randomly occurring red dots should not be evaluated in the interpretation of results.

Negative Result

A negative result is indicated by a horizontal line (Minus Sign “-”) in the Result Window. A negative result means that no hCG was detected, or that the levels of hCG in the specimen were below the detection limit of the assay.

Weak positive results may occur with hCG levels below 25mIU/ml. It is good laboratory practice to resample and retest these weak positive specimens after an additional 48-72 hours. The use of controls at hCG levels near the assay sensitivity may guide in the interpretation of weak positive results.

The reaction area may, on occasion, exhibit outlines. An outline can be described as a colourless area which surrounds all or part of the Patient Bar. If outlines are present, an impression of the Patient Bar may be seen. However, in the absence of hCG (negative specimens) this impression is comparable to the background and should be interpreted as a negative result.

Specimens with high levels of hCG may yield colour on the Patient Bar as quickly as 1 minute after specimen addition. Specimens with levels of hCG at or above the level of sensitivity of the assay remain positive with time. Specimens with levels of hCG below the sensitivity of the assay may yield some colour on the Patient Bar with time however the test should be read 5 minutes after addition of sample.

EXTERNAL QUALITY CONTROL

Good laboratory practice recommends the use of control materials to ensure proper kit performance. Each laboratory should refer to guidelines established internally and by local, national or other accrediting organisations.

Due to variation in analyte composition and/or matrices, external quality control materials and proficiency survey samples may not elicit identical results across all hCG assays. Each laboratory needs to determine the suitability of each control material for specific immunoassays and validate the material prior to use.

LIMITATIONS OF THE TEST

- Positive results from very early pregnancy may later prove negative due to natural termination of the pregnancy. This is estimated to occur in 31% of all conceptions⁶. It is recommended when using a sensitive pregnancy test such as Alere™ TestPack hCG Urine that weak positive results be retested with a first morning urine specimen taken 48-72 hours later.
- A negative result may be obtained if the urine specimen tested is too dilute.
- If a negative result is obtained and pregnancy is still suspected, the patient should be retested 48-72 hours later using a first morning urine specimen.
- Abnormal pregnancies (e.g. ectopic) may produce lower concentrations of hCG than expected for a given gestational age. Abnormal pregnancy cannot be distinguished from normal pregnancy by hCG levels alone^{7,8}.
- hCG remains elevated for a time after pregnancy⁹. Pregnancy tests carried out fewer than 3 weeks after giving birth, or 9 weeks after natural loss or termination, may need further evaluation.
- A number of conditions other than pregnancy can cause elevated levels of urinary hCG, e.g. menopause, ovarian cysts, trophoblastic disease, and certain non-trophoblastic neoplasms¹⁰.
- Occasionally, specimens containing <25mIU/ml hCG may test positive.
- Drugs containing hCG may interfere with Alere™ TestPack hCG Urine and produce misleading results.

9. False positive and false negative pregnancy tests may be observed in patients with abnormal bladder or kidney function e.g. enterocystoplasties and renal failure.
10. Inconsistent results may be obtained if the urine specimen contains an excessive amount of bacteria.
11. If the test result is not consistent with clinical evidence, further evaluation may be required.

EXPECTED RESULTS

Urine specimens from pre-menopausal females generally contain <5mIU/ml hCG; levels are generally <10mIU/ml in healthy males and post-menopausal females². On the first day of the first missed period, the levels of maternal hCG are normally 50-250mIU/ml³.

PERFORMANCE CHARACTERISTICS

Sensitivity

Alere™ TestPack hCG Urine can detect hCG in urine at concentrations of 25mIU/ml or greater. This sensitivity has been determined against the 4th International hCG Standard (WHO)⁴. Specimens containing less than 5mIU/ml should give negative results.

Prozone effect

Alere™ TestPack hCG Urine has been shown to produce positive results with specimens containing up to and including 1,000,000mIU/ml hCG, which is higher than the maximum level expected during typical pregnancy.

Specificity

Alere™ TestPack hCG Urine has been evaluated for its cross-reactivity with a variety of substances, including other hormones found in urine. No cross-reactivity was detected when the following substances were added to both "positive" (containing 25mIU/ml hCG) and "negative" urine specimens: LH (1000mIU/ml), FSH (1000mIU/ml), TSH (1000μIU/ml).

INTERFERING SUBSTANCES

The following potentially interfering substances were added to hCG negative and hCG positive (containing 25 mIU/ml hCG) specimens. None of the substances interfered with the assay at the concentrations tested.

Acetaminophen	(20 mg/dL)	Oxalic Acid	(60 mg/dL)
Acetoacetic Acid	(2000 mg/dL)	Oxytetracycline	(30 mg/dL)
Acetone	(1000 mg/dL)	Phenylpropanolamine	(4000 mg/dL)
Albumin (human serum)	(1200 mg/dL)	5β-Pregnane-3α, 20α-diol Glucuronide	(100 µg/dL)
Acetosalicylic Acid	(20 mg/dL)	Riboflavin	(2 mg/dL)
Ampicillin	(20 mg/dL)	Salicylic Acid	(20 mg/dL)
Ascorbic Acid	(200 mg/dL)	Sodium Carbonate	(800 mg/dL)
Atropine	(20 mg/dL)	Sodium Chloride	(6800 mg/dL)
Biotin	(25 µg/dL)	Tetracycline	(30 mg/dL)
Bilirubin	(1 mg/dL)	Urea	(2000 mg/dL)
Caffeine	(20 mg/dL)	Uric Acid	(100 mg/dL)
Creatinine	(200 mg/dL)		
Dextromethorphan	(20 mg/dL)	In addition no pH effect in range pH 4.5 to pH 8.5.	
Diphenhydramine	(20 mg/dL)		
EDTA	(40 mg/dL)		
Ephedrine	(20 mg/dL)		
Ethanol	(1%)		
Estrone β-D Glucuronide	(100 µg/dL)		
Gentisic Acid	(20 mg/dL)		
Glucose	(10000 mg/dL)		
Haemoglobin	(360 mg/dL)		
Hydroxybutyric Acid	(100 mg/dL)		
Human Serum Proteins	(2000 mg/dL)		
Ibuprofen	(40 mg/dL)		
Nicotine	(20 µg/dL)		

ACCURACY

In a study, 300 urine specimens collected from women for the purpose of pregnancy testing, were evaluated with the **Alere™ TestPack hCG Urine** and with the Clearview HCG Urine test.

Of the 300 urine specimens evaluated, 131 specimens tested positive by both methods and 169 specimens tested negative by both methods. Thus an agreement of 100% (300/300) was determined for these specimens. Therefore the diagnostic sensitivity and diagnostic specificity of **Alere™ TestPack hCG Urine** and Clearview HCG were found to be comparable in this study. As with all diagnostic tests, there is the potential for false positive or false negative results to occur with this test.

The study results are summarised below:

TestPack hCG Urine	Clearview HCG Urine	
	+	-
+	131	0
-	0	169

ADVICE LINE

For further information, please contact your distributor, or call Alere

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The Alere Logo and Alere are trademarks of the Alere group of companies.
© 2011 Alere. All rights reserved.
Alere Medical Co., Ltd.
Alere Ltd. Pepper Road, Hazel Grove,