Gebrauchsanweisung
Instructions for use
Instructions d'utilisation
Instrucciones de uso
Istruzioni per l'uso
Instrukcja obsługi

NADAL[®]

C.difficile Toxin A&B Test

(test cassette)

REF 582003

Version 1.0

2011-01-11

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14. Referenzen

- 1) Lierly, D. M., H.C. Krivan, and T.D. Wilkins. 1988. Clostridium difficile: Its disease and toxins. Clin. Microbiol. Rev. 1:1-18
- 2) Mulligan, M.E., L.R. Peterson, R.Y.Y. Kwok, C.R. Clabots, and D.N. Gerding. 1988. Immunoblots and plasmid fingerprints compared with serotyping and polyacrylamide gel electrophoresis for typing Clostridium difficile. J. Clin. Microbiol. 26:41-31
- 3) Gilligan, P.H., L.R. McCarthy, and V.M. Genta. 1981. Relative frequency of Clostridium difficile in patients with diarrheal disease. J. Clin. Microbiol. 26:41-31
- 4) George, W.L. 1989. Antimicrobial agent-associated diarrhea and colitis, p661-678. In S.M. Finegold and W.L. George (ed.), Anaerobic infections in humans. Academic Press, Inc. New York.
- 5) Sullivan, N.M., Pellet S., and T.D. Wilkins. 1982. Purification and characterization of toxins A and B of Clostridium difficile. Infect. Immun. 35:1032-1040
- 6) Lierly, D.M., D.E. Lockwood, S.H. Richardson, and T.D. Wilkins. 1982. Biological activities of toxins A and B of Clostridium difficile. Infect. Immun. 35:1147-1150

Ref.: 2011-01-11 SW

1. Intended Use

The NADAL[®] Clostridium difficile Toxin A&B Test is an immunochromatographic rapid test for the qualitative antigen detection of *Clostridium difficile* Toxin A and Toxin B in human feces.

2. Introduction and Clinical Significances

Clostridium difficile is a major cause of antibiotic-associated diarrhea and pseudomembranous colitis. It is one of the most commonly detected pathogens and an important cause of nosocomial infections in hospitals and nursing homes.

The organism has been isolated from diverse natural habitats, including soils, hay, sand, dung from various large mammals (cows, donkeys and horses), and from dog, cat, rodent and human feces. *Cl. difficile* produces at least three potential virulence factors from which Toxin A and Toxin B are thought to be the most important in the pathogenesis of *Cl. difficile* associated diseases.

Toxin A is an enterotoxin which seems to interfere with the cytoskeleton of the intestinal epithelial cells, rendering them non functional while Toxin B is a cytotoxin that induces strong cytopathic effects in tissue cultures cell lines.

Since not all strains of *Clostridium difficile* produce toxins and approx. 2% of healthy adults as well as up to 50% of children younger than 2 years can be colonized with *Clostridium difficile*, the detection of the toxins (Toxin A and Toxin B) in stool samples of patients with diarrhea is more significant than culturing the bacteria.

3. Principle of the Test

The test device consists of a plastic housing containing two different sticks for the detection of *Cl. difficile* Toxin A or *Cl. difficile* Toxin B. A coloured anti *Cl. difficile* Toxin A or *Cl. difficile* Toxin B monoclonal antibody colloidal gold conjugate is placed at the left end of the membrane.

After collection in a tube containing the extraction solution, the feces sample is dissolved and few drops of this extract are added into each well of the reaction device.

As the test sample flows through the absorbent device, the labelled antibody-dye conjugate binds to the Toxin A or Toxin B antigen (when present in the sample), forming an antibody antigen complex.

This complex binds to the polyclonal anti-Toxin A or anti-Toxin B antibody in the positive reaction zone, producing a rose-pink coloured band.

In the absence of Toxin A or Toxin B, there is no line in the positive reaction zone. The reaction mixture continues flowing through the absorbent device, past the positive reaction zone and control zone.

Unbound conjugate binds to the reagent in the control zone producing a rose-pink coloured band demonstrating that the reagents are functioning correctly.

4. Reagents and Materials Supplied

Each NADAL[®] Clostridium difficile Toxin A&B Testkit contains everything needed to perform 10 tests:

- 10 NADAL[®] Clostridium difficile Toxin A&B test cassettes
- 10 single use plastic pipettes
- 10 extraction tubes (filled with 2 mL of extraction solution)
- 10 sample applicators
- 1 instructions leaflet

5. Additional Required Materials

- Timer
- Stool collector

6. Storage & Stability

All NADAL[®] Clostridium difficile Toxin A&B Test kit components should be stored in the sealed pouch at room temperature (between +2°C and +30°C). Do not freeze the test kit. The NADAL[®] Clostridium difficile Toxin A&B Test kit is stable until the expiry date stated on the package label.

7. Warnings and Precautions

- This test is designed for *in vitro* diagnostic use and professional use only.
- Read carefully instructions leaflet before using this test.
- Do not use beyond the expiry date stated on the package label.
- Do not use a test from a damaged protective wrapper.
- All reagents and materials coming in contact with potential infectious specimens must be treated with appropriate disinfectants or autoclaved at 121°C for at least one hour.
- Wear protective clothing such as laboratory coats and disposable gloves while assaying samples.
- Do not eat, drink or smoke in the area where specimens and kit reagents are handled.
- Avoid any contact between hands and eyes or nose during specimen collection and testing.

8. Specimen Collection and Preparation

Preliminary notes for specimen collection and storage

Stool specimen should be collected as soon as possible after onset of symptoms.

Diluted samples may be stored at +2°C to +8°C in refrigerator for 3 days without interference with assay performance.

For long term storage of undiluted specimens, storage at -20°C or colder is recommended. Repeated freezing and thawing of samples is not recommended and may cause erroneous results.

Caution!

Do not collect specimens in containers having media, preservatives, animal serum or detergents, as any of these may interfere with the test.


Procedure of specimen collection

- 1) Label the plastic tube containing the extraction solution with patient's name or control number.
- 2) Open the extraction tube.
- 3) Transfer stool sample into the extraction tube, using sample applicator:
 - volume of a pea, in case of solid stool
- 4) Close the extraction tube by tighten the red cap of the tube.
- 5) Mix the stool sample and the diluent by shaking well until the sample is dissolved.
- 6) Let the tube stand long enough for the large particles to settle to the bottom of the tube. Alternatively, centrifuge the tube at 500-1,000 RPM for 1 minute.

9. Procedure of the Test

- 1) Bring all reagents at room temperature (15-30°C).
- 2) Remove the test device from the pouch by tearing along the split.
- 3) Open the plastic tube containing the extracted sample.
- 4) Using the plastic pipette, add 6 drops (200 µl) of the extracted solution into each of the sample wells for TOX-A and TOX-B.
- 5) Read the test results 15 minutes after addition of the sample to the device.

10. Interpretation of the Results

Negative

Only one coloured band appears in the control zone. No band is visible in the test zone.


Positive

In addition to the control band, a clearly distinguishable band also appears in the test zone.

Invalid

If there is no distinct coloured band visible in the control zone, the test is inconclusive. Repeat the test.

11. Quality Control

A clear background in the observation window is considered an internal negative control. However, when the faecal sample is tested, the background may appear slightly reddish due to the original colour of the sample. This is acceptable as long as it does not interfere with the interpretation of test result. The test is invalid if the background fails to clear and obscures the reading of the result.

12. Limitations

- The NADAL[®] Clostridium difficile Toxin A&B Test is specifically designed to detect Toxin A or Toxin B antigen in the stool samples.
- As for any *in vitro* diagnostic procedure, the physician should confirm the test results with other clinical methods.
- A negative result does not generally exclude a *Cl.difficile* infection. It can be caused by proteolytic digestion of the toxins due to improper specimen storage. If a reasonable suspicion of an infection exists, another stool specimen should be investigated.
- A positive result does not exclude the presence of other pathogens.
- The presence of blood in the feces samples in significant quantity may lead to false positive results in limited cases.
- Test and control lines colours may slightly change depending on the stool sample aspect. For example dark green lines (instead of pink lines) have been reported when assaying greenish or darkish stool samples. This stool coloration appears in case of treatment of iron deficiency with ferrous fumarate. The test result should be interpreted as usual, i.e. two lines for a positive result and one line for a negative result.

13. Performance Characteristics

Cl.difficile Toxin A

Analytical sensitivity

The performance of the NADAL[®] Clostridium difficile Toxin A&B Test has been assayed, using a range diluted solutions prepared from a commercially available purified toxin A antigen. The NADAL[®] Clostridium difficile Toxin A&B Test is able to detect a concentration of 4 ng/mL.

Diagnostic sensitivity and specificity

A study was performed on 63 feces samples, using the NADAL[®] Clostridium difficile Toxin A&B Test test in comparison with the cytotoxicity and cell culture methods.

The results are summarized in the table 1:

		NADAL [®] Clostridium difficile Toxin A&B		Total
		+	-	
Cytotoxicity	+	13	2	15
	-	4	44	48
Total		17	46	63

All samples showed identical results using the cytotoxicity and the cell culture methods.

From the above table, the sensitivity of NADAL[®] Clostridium difficile Toxin A&B Test test is 86.7% (13/15), compared to the cytotoxicity method and the specificity is 91.7% (44/48), compared to the cytotoxicity method.

Cross reaction

NADAL[®] Clostridium difficile Toxin A&B Test showed consistently negative results up to 500 ng/mL Toxin B.

Cl.difficile Toxin B

Analytical sensitivity

The sensitivity of the NADAL[®] Clostridium difficile Toxin A&B Test has been evaluated using a range of diluted solutions prepared from a commercially available purified *C. difficile* Toxin B antigen. Under these conditions, the detection limit of the NADAL[®] Clostridium difficile Toxin A&B Test has been found to be 5 ng/mL.

Specificity

A purified *Cl. difficile* Toxin A antigen was used to determine the specificity of the test. NADAL[®] Clostridium difficile Toxin A&B Test showed consistently negative results up to 500 ng/mL Toxin A.

14. References

- 1) Lierly, D, M, H.C. Krivan, and T.D. Wilkins. 1988. Clostridium difficile: Its disease and toxins. Clin. Microbiol. Rev. 1:1-18
- 2) Mulligan, M.E., L.R. Peterson, R.Y.Y. Kwok, C.R. Clabots, and D.N. Gerding. 1988. Immunoblots and plasmid fingerprints compared with serotyping and polyacrylamide gel electrophoresis for typing Clostridium difficile. J. Clin. Microbiol. 26:41-31
- 3) Gilligan, P.H., L.R. McCarthy, and V.M. Genta. 1981. Relative frequency of Clostridium difficile in patients with diarrheal disease. J. Clin. Microbiol. 26:41-31
- 4) George, W.L. 1989. Antimicrobial agent-associated diarrhea and colitis, p661-678. In S.M: Finegold and W.L. George (ed.), Anaerobic infections in humans. Academic Press, Inc. New York.
- 5) Sullivan, N.M., Pellet S., and T.D. Wilkins. 1982. Purification and characterization of toxins A and B of Clostridium difficile. Infect. Immun. 35:1032-1040
- 6) Lyerly, D.M., D.E. Lockwood, S.H. Richardson, and T.D. Wilkins. 1982. Biological activities of toxins A and B of Clostridium difficile. Infect. Immun. 35:1147-1150

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14. Références

- Lierly, D. M., H.C. Krivan, and T.D. Wilkins. 1988. Clostridium difficile: Its disease and toxins. Clin. Microbiol. Rev. 1:1-18
- Mulligan, M.E., L.R. Peterson, R.Y.Y. Kwok, C.R. Clabots, and D.N. Gerding. 1988. Immunoblots and plasmid fingerprints compared with serotyping and polyacrylamide gel electrophoresis for typing Clostridium difficile. J. Clin. Microbiol. 26:41-31
- Gilligan, P.H., L.R. McCarthy, and V.M. Genta. 1981. Relative frequency of Clostridium difficile in patients with diarrheal disease. J. Clin. Microbiol. 26:41-31
- George, W.L. 1989. Antimicrobial agent-associated diarrhea and colitis, p661-678. In S.M: Finegold and W.L. George (ed.), Anaerobic infections in humans. Academic Press, Inc. New York.
- Sullivan, N.M., Pellet S., and T.D. Wilkins. 1982. Purification and characterization of toxins A and B of Clostridium difficile. Infect. Immun. 35:1032-1040
- Lyerly, D.M., D.E. Lockwood, S.H. Richardson, and T.D. Wilkins. 1982. Biological activities of toxins A and B of Clostridium difficile. Infect. Immun. 35:1147-1150

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13. Referencias

1. Lierly, D.M, H.C. Krivan, and T.D. Wilkins. 1988. Clostridium difficile : its disease and toxins. Clin. Microbiol. Rev. 1:1-18.
2. Mulligan, M.E., L.R. Peterson, R.Y.Y. Kwok, C.R. Clabots, and D.N. Gerding. 1988. Immunoblots and plasmid fingerprints compared with serotyping and polyacrylamide gel electrophoresis for typing Clostridium difficile. J. Clin. Microbiol. 26 : 41-31.

3. Gilligan, P.H., L.R. McCarthy, and V.M. Genta. 1981. Relative frequency of Clostridium difficile in patients with diarrheal disease. J. Clin. Microbiol. 14 : 26-31.
4. George, W.L. 1989. Antimicrobial agent-associated diarrhea and colitis, p661-678. In S.M. Finegold and W.L. George (ed.), Anaerobic infections in humans. Academic Press, Inc. New York.
5. Sullivan, N.M., Pellet S., and T.D. Wilkins. 1982. Purification and characterization of toxins A and B of Clostridium difficile. Infect. Immun. 35 : 1032-1040.
6. Lierly, D.M., D.E. Lockwood, S.H. Richardson, and T.D. Wilkins. 1982. Biological activities of toxins A and B of Clostridium difficile. Infect. Immun. 35 : 1147-1150.

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14. Bibliografia

- 7) Lierly, D. M., H.C. Krivan, and T.D. Wilkins. 1988. Clostridium difficile: Its disease and toxins. Clin. Microbiol. Rev. 1:1-18
- 8) Mulligan, M.E., L.R. Peterson, R.Y.Y. Kwok, C.R. Clabots, and D.N. Gerding. 1988. Immunoblots and plasmid fingerprints compared with serotyping and polyacrylamide gel electrophoresis for typing Clostridium difficile. J. Clin. Microbiol. 26:41-31
- 9) Gilligan, P.H., L.R. McCarthy, and V.M. Genta. 1981. Relative frequency of Clostridium difficile in patients with diarrheal disease. J. Clin. Microbiol. 26:41-31
- 10) George, W.L. 1989. Antimicrobial agent-associated diarrhea and colitis, p661-678. In S:M: Finegold and W.L. George (ed.), Anaerobic infections in humans. Academic Press, Inc. New York.
- 11) Sullivan, N.M., Pellet S., and T.D. Wilkins. 1982. Purification and characterization of toxins A and B of Clostridium difficile. Infect. Immun. 35:1032-1040
- 12) Lierly, D.M., D.E. Lockwood, S.H. Richardson, and T.D. Wilkins. 1982. Biological activities of toxins A and B of Clostridium difficile. Infect. Immun. 35:1147-1150

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

14. Bibliografia

1. Lierly, D, M, H.C. Krivan, and T.D. Wilkins. 1988. Clostridium difficile: Its disease and toxins. Clin. Microbiol. Rev. 1:1-18
2. Mulligan, M.E., L.R. Peterson, R.Y.Y. Kwok, C.R. Clabots, and D.N. Gerding. 1988. Immunoblots and plasmid fingerprints compared with serotyping and polyacrylamide gel electrophoresis for typing Clostridium difficile. J. Clin. Microbiol. 26:41-31
3. Gilligan, P.H., L.R. McCarthy, and V.M. Genta. 1981. Relative frequency of Clostridium difficile in patients with diarrheal disease. J. Clin. Microbiol. 26:41-31
4. George, W.L. 1989. Antimicrobial agent-associated diarrhea and colitis, p661-678. In S:M: Finegold and W.L. George (ed.), Anaerobic infections in humans. Academic Press, Inc. New York.
5. Sullivan, N.M., Pellet S., and T.D. Wilkins. 1982. Purification and characterization of toxins A and B of Clostridium difficile. Infect. Immun. 35:1032-1040
6. Lyerly, D.M., D.E. Lockwood, S.H. Richardson, and T.D. Wilkins. 1982. Biological activities of toxins A and B of Clostridium difficile. Infect. Immun. 35:1147-1150

Symbol	English	Deutsch	Français	Nederlands	Español	Italiano
	Consult instructions for use	Gebrauchsanweisung beachten	Consulter les instructions d'utilisation	Gebruiksaanwijzing	Consulte las instrucciones de uso	Consultare le istruzioni per l'uso
	European Conformity	CE-Konformitätskennzeichnung	Conformité aux normes européennes	Conform Europese richtlijnen	Conformidad europea	Conformità europea
	In vitro diagnostic device	In-vitro-Diagnostikum	Usage Diagnostic in vitro	In vitro diagnostisch gebruik	Para uso Diagnóstico in vitro	Per uso Diagnostica in vitro
	Catalogue number	Katalog-Nr.	Numéro de catalogue	Catalogus nummer	Número de catálogo	Numero di Catalogo
	Lot. No. / Batch code	Chargen-Nr.	Numéro de lot	Lot nummer	Número de lote	Numero di lotto
	Contains sufficient for <n> tests/	Ausreichend für "n" Ansätze	Contenu suffisant pour "n" tests	Geschikt voor <n> tests	Contenido suficiente para <n> ensayos	Contenuto sufficiente per "n" saggi
	Storage Temperature	Lagerungs-temperatur	Température de conservation	Opslag temperatuur	Temperatura de conservación	Temperatura di conservazione
	Expiration Date	Mindesthaltbarkeitsdatum	Date limite d'utilisation	Houdbaarheids datum	Fecha de caducidad	Data di scadenza
	nal von minden GmbH, Carl-Peschken-Straße 9, 47441 Moers, Germany	nal von minden GmbH, Carl-Peschken-Straße 9, 47441 Moers, Germany	nal von minden GmbH, Carl-Peschken-Straße 9, 47441 Moers, Germany	nal von minden GmbH, Carl-Peschken-Straße 9, 47441 Moers, Germany	nal von minden GmbH, Carl-Peschken-Straße 9, 47441 Moers, Germany	nal von minden GmbH, Carl-Peschken-Straße 9, 47441 Moers, Germany
	Content	Inhalt	Conditionnement	Inhoud	Contenido	Contenuto

Symbol	Polski	Suomi	Portugues	Dansk	Svenska	Ελληνικά
	Patrz: ulotka informacyjna	Katso käyttöohjeet	Consulte as instruções de utilização	Se brugsanvisning	Se bruksanvisningen	Εγχειρίδιο χρήστη
	Znak zgodności CE	CE-merkitty	Conformidade com as normas europeias	Europeaisk overensstemmelse	Europeisk överensstämmelse	Ευρωπαϊκή Συμμόρφωση
	Tylko do użytku in-vitro	In vitro-diagnostiikka	Diagnóstico in vitro	In vitro diagnostik	In vitro diagnostik	in vitro διαγνωστικό
	Numer katalogowy	Luettelonumero	Catálogo n.º	Katalognummer	Katalognummer	Αριθμός καταλόγου
	Numer serii	Eränumero	No do lote	Lot nummer	Batchnummer	Αριθμός Παρτίδος
	Wystarszające na "n" powtórzeń	Sisältää tarvikkeet "n" testiin		Indeholder tilstrækkeligt til "n" test	Innehåller tillräckligt till "n" tester	Περιεχόμενο επαρκές για «n» εξετάσεις
	Temperatura przechowywania	Säilytyslämpötila	Temperatura de conservação	Opbevarings-temperatur	Förvaringstemperatur	Θερμοκρασία αποθήκευσης
	Data ważności	Viimeinen käyttöpäivä	Prazo de validade	Udløbsdato	Utgångsdatum	Ημερομηνία λήξης
	nal von minden GmbH, Carl-Peschken-Straße 9, 47441 Moers, Germany	nal von minden GmbH, Carl-Peschken-Straße 9, 47441 Moers, Germany	nal von minden GmbH, Carl-Peschken-Straße 9, 47441 Moers, Germany	nal von minden GmbH, Carl-Peschken-Straße 9, 47441 Moers, Germany	nal von minden GmbH, Carl-Peschken-Straße 9, 47441 Moers, Germany	nal von minden GmbH, Carl-Peschken-Straße 9, 47441 Moers, Germany
	Zawartość	Sisältö	Conteúdo	Indhold	Innehåll	Περιεχόμενο

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