HIV Combo

Key to symbols used / Schlüssel für die verwendeten Symbole / Clave de los símbolos utilizados / Légende des symboles utilisés / Legenda dei simboli utilizzati / Descrição dos símbolos utilizados

 **2°C - 30°C** Store at 2-30°C / Bei 2-30°C lagern / Guardar a temperaturas entre 2 y 30°C / Conserver entre 2-30°C / Conservare a 2-30°C / Armazenar entre 2-30°C

EDTA CAPILLARY TUBES

EDTA Capillary Tubes / EDTA-Kapillarröhrchen / Tubos capilares EDTA / Tubes capillaires avec EDTA / Tubi capillari con EDTA / Tubos Capilares EDTA

CHASE BUFFER

Chase Buffer / Chasepuffer / Buffer de detección / Tampon de migration / Tampone chase (tampone di spinta) / Tampão de detecção

Alere Ltd.
Pepper Road, Hazel Grove,
Stockport, SK7 5BW, UK

Alere Medical Co., Ltd.
357 Matsuhidai, Matsudo-shi,
Chiba, 270-2214 Japan
Tel +81 47 311 5750

www.alere.com

0088

Alere™ HIV Combo

EN

This package insert must be read carefully prior to use. Package insert instructions must be followed accordingly. Reliability of assay results cannot be guaranteed if there are deviations from the instructions in this package insert.

NAME AND INTENDED USE

Alere™ HIV Combo is an *In Vitro*, visually read, qualitative immunoassay for the simultaneous detection of free non immunocomplexed HIV-1 p24 antigen (Ag) and antibodies (Ab) to HIV-1 and HIV-2 in human blood. The test specimen can be serum, plasma, fingerstick or venous whole blood. The test is intended as a rapid test to detect free non immunocomplexed HIV-1 p24 antigen and antibodies to HIV-1/HIV-2 from infected individuals.

SUMMARY AND EXPLANATION OF THE TEST

AIDS (Acquired Immunodeficiency Syndrome) is characterized by changes in the population of T-cell lymphocytes. In an infected individual, the virus causes depletion of helper T-cells, which leaves the person susceptible to opportunistic infections and some malignancies. The virus that causes AIDS exists as two related types known as HIV-1 and HIV-2. The presence of HIV first elicits the secretion of p24 antigen^{1,2} followed by the production of specific antibodies to either HIV-1 or HIV-2.^{3,4,5}

BIOLOGICAL PRINCIPLES OF THE PROCEDURE

Alere™ HIV Combo is an immunochemical test for the qualitative detection of HIV-1 p24 antigen and antibodies to HIV-1 and HIV-2. Specimen is added to the sample pad. The specimen mixes with biotinylated anti-p24 antibodies and selenium colloid - conjugates coated with recombinant HIV-1, HIV-2 and HIV-1 group O antigens, synthetic HIV-2 peptide and anti p24 mouse monoclonal antibody. This mixture continues to migrate through the solid phase to the immobilized recombinant HIV-1/HIV-1 group O antigens and synthetic HIV-1/HIV-2 peptides at the Antibody (Ab) window, immobilized avidin at the Antigen (Ag) window.

If antibodies to HIV-1 and/or HIV-2 are present in the specimen, the antibodies bind to the selenium colloid- conjugates coated with recombinant HIV-1, HIV-2 and HIV-1 group O antigens and synthetic HIV-2 peptide and to the immobilized recombinant HIV-1/HIV-1 group O antigens and synthetic HIV-1/HIV-2 peptides, forming one red bar at the Ab window site. If antibodies to HIV-1 and HIV-2 are absent, the selenium colloid- conjugates flow past the Ab window and no red bar is formed at the Ab window site.

If free non immunocomplexed HIV-1 p24 antigen is present in the specimen, the antigen binds to the biotinylated anti-p24 antibodies and the selenium colloid-conjugate coated with anti p24 mouse monoclonal antibody and it binds to an immobilized avidin forming a red bar at the Ag window site. If HIV-1 p24 antigen is not present, both the biotinylated anti-p24 antibodies and selenium colloid-conjugate flow past the Ag window and no red bar is formed at the Ag window site.

To ensure assay validity, a procedural control bar is incorporated in the assay device at the Control window.

CONTENTS

Alere™ HIV Combo 20 Test (7D2846) or 100 Test (7D2847)

- Alere™ HIV Combo Test Card, 2 or 10 cards (10 tests/card) coated with HIV-1/2 recombinant antigen, synthetic peptides, anti-p24 antibodies and avidin.

ACCESSORIES (required but not provided)

For testing Whole Blood samples

CHASE BUFFER | 1 Bottle (2.5 mL) Chase Buffer (7D2243) prepared in phosphate buffer.

Preservatives: Antimicrobial Agents.

Whole Blood (fingerstick assay)

EDTA CAPILLARY TUBES | (7D2227)

MICROSAFE® Tube (SAFE-TEC® Clinical Products, YY1050 (2000 tubes), YY1050100 (100 tubes))

Lancet

WARNINGS AND PRECAUTIONS

For *In Vitro* Diagnostic Use.

For Professional use only.

Safety data sheet available for professional user on request.

Patients with elevated triglycerides levels may test false reactive with the **Alere™ HIV Combo**.

CAUTION:

Appropriate biosafety practices^{7,8} should be used when handling specimens and reagents. These precautions include, but are not limited to the following:

- Wear gloves.
- Do not pipette by mouth.
- Do not eat, drink, smoke, apply cosmetics, or handle contact lenses in areas where these materials are handled.
- Clean and disinfect all spills of specimens or reagents using suitable disinfectant, such as 0.5% sodium hypochlorite, or other suitable disinfectant.⁹
- Decontaminate and dispose of all specimens, reagents, and other potentially contaminated materials in accordance with local regulations.^{10,11}

STORAGE

Alere™ HIV Combo Test Cards and Chase Buffer must be stored at 2-30°C until expiration date.

- Kit components are stable until the expiration date when handled and stored as directed. Do not use kit components beyond expiration date.
- Immediately reseal all unused tests in the foil pouch containing the desiccant by pressing seal from end to end to close.
- Do not use devices that have become wet or if the packaging has become damaged.

SPECIMEN COLLECTION

Serum, Plasma, and Whole Blood Collection by Venipuncture.

- Human serum, plasma, and whole blood collected by venipuncture should be collected aseptically in such a way as to avoid hemolysis.
- Separate the serum from the clot or plasma from the packed cells as soon as possible to avoid any hemolysis.
- Specimens showing particulate matter or turbidity should be centrifuged before testing in order to avoid providing erroneous results. The specimen must be centrifuged at 10,000g for 5 min at room temperature. Carefully remove the 50 µL test sample from the supernatant. If a lipid layer is formed on the surface of the liquid, ensure that the sample is taken from the clear liquid below that layer.

NOTE: For whole blood and plasma specimens, EDTA collection tubes must be used.

Whole Blood Collection by Fingerstick¹²

EDTA CAPILLARY TUBES or **MICROSAFE® Tube** must be used.

Before collecting a fingerstick specimen, place an EDTA capillary tube on a clean dry surface.

1. Choose the fingertip of the middle, ring, or index finger (whichever is the least callused). Warm the hand as needed with a warm, moist towel or warm water to increase blood flow.
2. Clean fingertip with alcohol; allow to air dry.
3. Position the hand palm-side up. Place the lancet off-center on the fingertip. Firmly press the lancet against the finger and puncture the skin. Dispose of the lancet in an appropriate biohazard sharps container.
4. Wipe away the first drop of blood with a sterile gauze pad.
5. Hold the finger lower than the elbow and apply gentle, intermittent pressure to the base of the punctured finger several times. Touch the tip of the EDTA Capillary Tube to the drop of blood¹³. Avoid air bubbles.
6. Fill the tube with blood between the 2 marked lines. (50 µL).

* If **MICROSAFE® Tube** will be used, refer to **MICROSAFE® Tube** package insert for additional information.

SPECIMEN STORAGE

- Serum and plasma specimens should be stored at 2-8°C if the test is to be run within 7 days of collection. If testing is delayed more than 7 days, the specimen should be frozen (-20°C or colder).

• Avoid repeated freeze/thaw cycles.

• Whole blood collected by venipuncture should be stored at 2-8°C if the test is to be run within 7 days of collection. Do not freeze whole blood specimens.

Mix the specimen well by gentle inversion of the tube immediately before testing.

• Whole blood collected by fingerstick should be tested immediately.

TEST PROCEDURE

The desired number of test units from the 10-test card can be removed by bending and tearing at the perforation.

NOTE: Removal of the test units should start from the right side of the test card to preserve the lot number which appears on the left side of the test card.

• The assay should be initiated immediately after removing the protective foil cover from each test.

• If the test strip dries during reading and it becomes difficult to see the bars, the test should be repeated.

• If the sample does not flow or shows abnormal flow, such as stopping in the middle of the window, the specimen should be centrifuged and the test repeated.

1. Remove the protective foil cover from each test.

2. For serum or plasma samples:

- a. Apply 50 µL of sample (precision pipette) to the sample pad (marked by the arrow symbol).
- b. Wait a minimum of 20 minutes from addition of the sample (up to 40 minutes maximum) and read result.

3. For whole blood (venipuncture) samples and for whole blood (fingerstick) samples using a **MICROSAFE® Tube**:

- a. Apply 50 µL of sample (precision pipette or **MICROSAFE® Tube**) to the sample pad (marked by the arrow symbol).
- b. Wait one minute to allow the sample to be absorbed, and then apply one drop of Chase Buffer to the sample pad.
- c. Wait a minimum of 20 minutes from addition of the sample (up to 40 minutes maximum) and read result.

4. For whole blood (fingerstick) samples using an EDTA Capillary Tube:

- a. Place the capillary tube containing the blood sample onto the sample pad (marked by the arrow symbol).
- b. Wait until all the blood is transferred from the capillary tube to the sample pad and then immediately apply one drop of Chase Buffer to the sample pad. Caution: do not lift the capillary from the sample pad before all the blood has been transferred – a bubble may form which will prevent the complete transfer of sample and invalidate the test. It may take more than one minute for full transfer of the sample.
- c. Wait a minimum of 20 minutes from addition of the sample (up to 40 minutes maximum) and read result.


QUALITY CONTROL

To ensure assay validity, a procedural control is incorporated in the device and is labeled "C". Any visible bar (even very faint) in the control window should be interpreted as a valid result. If the control bar does not appear by assay completion, the test result is invalid and the sample should be retested.

INTERPRETATION OF RESULTS

ANTIBODY POSITIVE (Two Bars - Control and Ab Bars)

Red bars appear in both the control window (labeled "C") and in the Ab window (labeled "AB") of the strip. Any visible red bar in the Ab window should be interpreted as positive. The presence of only an antigen response suggests that the infection is at an early stage. Follow up testing may be suggested in order to track the expected future detection of antibodies.

ANTIGEN (p24) POSITIVE (Two Bars - Control and Ag Bars)

Red bars appear in both the control window (labeled "C") and in the Ag window (labeled "AG") of the strip. Any visible red bar in the Ab and Ag windows should be interpreted as positive.

ANTIBODY POSITIVE AND ANTIGEN (p24) POSITIVE (Three Bars - Control, Ab and Ag Bars)

Red bars appear in the control window (labeled "C"), the Ab window (labeled "AB") and the Ag window (labeled "AG") of the strip. Any visible red bar in the Ab and Ag windows should be interpreted as positive.


NEGATIVE (One Bar - Control Bar)

One red bar appears in the control window of the strip (labeled "C"), and no red bar appears in the Ab and Ag windows of the strip (labeled "AG" and "AB").

INVALID (No Control Bar)

If there is no red bar in the control window of the strip, and even if a red bar appears in the Ab or Ag window of the strip, the result is invalid and should be repeated.

NOTES:

- The test result is positive even if the patient bars appear lighter or darker than the control bar.
- If an invalid result occurs repeatedly, or for technical assistance, contact your local distributor or call Technical Support.

LIMITATIONS OF THE PROCEDURE

- **Alere™ HIV Combo** is designed to simultaneously detect antibodies to HIV-1 and/or HIV-2 and free non immunocomplexed HIV-1 p24 antigen, in human serum, plasma and whole blood. Other body fluids or pooled specimens may not give accurate results and should not be used.
- The intensity of the Ab and Ag bars does not correlate to the titer of antibody and antigen in the specimen.
- No test provides absolute assurance that a specimen does not contain low levels of HIV-1 p24 antigen and/or antibodies to HIV-1 and HIV-2 such as those present at a very early stage of infection.
- A negative result for both antibodies to HIV and HIV-1 p24 antigen does not preclude the possibility of exposure to or infection with HIV-1 or HIV-2 viruses.
- Where clinical presentation or other data would suggest an inconsistent result then the patient should be tested by PCR immediately and/or retested for antibodies to HIV >21 days after the original testing.
- A positive result for antibodies to HIV with a negative result for HIV-1 p24 antigen does not preclude the possibility of acute infection.
- Positive results should be confirmed using another method and the results should be evaluated in light of the overall clinical evaluation before a diagnosis is made.
- Whole blood or plasma specimens containing anticoagulants other than EDTA have not been validated for use with the **Alere™ HIV Combo** and may give incorrect results.
- Infants born to HIV-infected mothers may carry maternal antibodies and will test antibody positive until eighteen months of age, which may not necessarily indicate the true infection status of the newborn. The use of HIV-1 p24 antigen testing to exclude infection in neonates (up to around eighteen months) is not recommended by CDC, because of poor sensitivity, especially in the presence of HIV antibody.
- Definitive diagnosis of HIV infection in early infancy requires other assays, including HIV nucleic acid test or viral culture.⁹
- Some known HIV-infected persons taking Antiretroviral medication have been shown to produce false negative results when tested by rapid kits.^{10,11,12}

PERFORMANCE CHARACTERISTICS

SENSITIVITY

Sensitivity was evaluated by testing confirmed HIV antibody positive specimens, commercial seroconversion panels and HIV-1 subtype virus panels.

1. HIV Antibody Positive Specimens:

Table I
HIV Antibody Positive Specimens

Types	Number of Specimens Tested	Positive by Alere™ HIV Combo	Sensitivity
HIV-1	422	422	100.0%
HIV-1 non B subtypes *	56	56	100.0%
HIV-1 group O	4	4	100.0%
HIV-2	100	100	100.0%
HIV**	100	100	100.0%
Total	682	682	100.0%

* Subtypes: A, C, D, F, G, H, J, K, and CRF01-AE, CRF02-AG, CRF03-AB, CRF05-DF, CRF09-A/U, CRF11-cpx

** Specimens collected in Europe

A total of 682 confirmed HIV positive specimens were tested (Table I). The antibody diagnostic sensitivity of **Alere™ HIV Combo** on this population of specimens was calculated to be 100%.

2. HIV Seroconversion Panels:

Table II
HIV Seroconversion Panels

Compared to	No. of seroconversion panels	Alere™ HIV Combo results compared to each product		
		Earlier detection (one bleed or more)	Equivalent detection (Same bleed)	Later detection (one bleed later)
CE marked 3 rd generation rapid test	31	28	3	0
CE marked 4 th generation rapid test	31	10	21	0
CE marked 4 th generation kit (CMIA)	20	0	14	6

The sensitivity of the **Alere™ HIV Combo** was evaluated with 31 seroconversion panels. The results were compared with the results of commercially available CE marked HIV 3rd generation rapid test kit and CE marked HIV 4th generation rapid test kit. Additionally, the results of 20 seroconversion panels were compared with the data of CE marked HIV 4th generation kit (Chemiluminescent microparticle immunoassay (CMIA)). **Alere™ HIV Combo** detected HIV infection earlier than 3rd generation rapid test kit and 4th generation rapid test kit. Detection of six panels was delayed by 1 bleed date from 4th generation kit (CMIA).

3. HIV-1 p24 Subtype Antigens:

Table III
HIV-1 p24 Subtype Antigens

Ag Subtype	Number of virus panels tested	Antigen Positive by Alere™ HIV Combo		Sensitivity
HIV-1 A	4	4		100.0%
HIV-1 B	8	8		100.0%
HIV-1 C	7	7		100.0%
HIV-1 D	5	5		100.0%
HIV-1 AE	10	10		100.0%
HIV-1 AG	2	2		100.0%
HIV-1 F	5	5		100.0%
HIV-1 G	2	2		100.0%
HIV-1 H	1	1		100.0%
HIV-1 Group O	4	4		100.0%
HIV-1 Total	48	48		100.0%
HIV-2	2	0		0.0%

50 HIV virus panels were prepared by diluting 50 different cultured HIV isolates including HIV-1 group M subtypes: A, B, C, D, F, G, H and CRF01-AE, CRF02-AG, HIV-1 group O and HIV-2 with HIV negative human serum. The **Alere™ HIV Combo** could detect all HIV-1 panels and could not detect HIV-2 panels (Table III).

4. Analytical Sensitivity of HIV-1 p24 antigen

The analytical sensitivity of the **Alere™ HIV Combo** was evaluated by testing WHO International Standard HIV-1 p24 Antigen (NIBS code 90/636). The results demonstrated that the test could detect a concentration of 2 IU/mL HIV-1 p24 antigen.

The antigen detection sensitivity of **Alere™ HIV Combo** with 29 seroconversion panels based on RNA copies/ml was 14.3% (2/14) in the range of 1.0x10⁴ to 9.9x10⁴ copies/ml, 70.8% (17/24) in the range of 1.0x10⁵ to 9.9x10⁵ copies/ml, and 100% (6/6) in the range of 1.0x10⁶ to 9.9x10⁶ copies/ml.

SPECIFICITY

A total of 2469 confirmed negative serum or plasma specimens were tested with the **Alere™ HIV Combo** and the specificity was determined for the Ab bar and for the Ag bar (Table IV). The specificity was 99.95% for the Ab bar and 99.76% for the Ag bar.

Table IV
Specificity of Alere™ HIV Combo
Population	Number of Tested Specimens	Alere™ HIV Combo			
Ab bar	Ag bar	Combined Ab and Ag bar			

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

BIBLIOGRAPHY

1. Pilcher C, Eron JJ, Galvin S, et al. (2004) Acute HIV revisited: new opportunities for treatment and prevention. *The Journal of Clinical Investigations* 113(7): 937-945.
2. Respass RA, Rayfield MA and Dondero TJ (2001) Laboratory testing and rapid HIV assays: applications for HIV surveillance in hard-to-reach populations. *AIDS* 15 Supplement 3: S49-S59.
3. Louie B, Wong E, Klausner JD, et al. (2008) Assessment of rapid tests for detection of Human Immunodeficiency Virus-specific antibodies in recently infected individuals. *Journal of Clinical Microbiology* 46(4): 1494-1497.
4. Lyons MS, Lindsell CJ, Hawkins DA, et al. (2008) Contributions to early HIV diagnosis among patients linked to care vary by testing venue. *BMC Public Health* 8:220
5. Gürler LG, Hauser PH, Eberle J, et al. A New Subtype of Human Immunodeficiency Virus Type 1 (MVP-5180) from Cameroon. *Journal of Virology*. 1994; 68(3):1581-1585.
6. Clinical and Laboratory Standards Institute. Procedures and Devices for the Collection of Diagnostic Capillary Blood Specimens; Approved Standard-Sixth Edition GP42-A6 Vol.28 No.25 September 2008
7. Clinical and Laboratory Standards Institute. Clinical Laboratory Waste Management; Approved Guideline -Third Edition. GP05-A3 Vol.31 No.3 January 2011
8. EPA Guide for Infection Waste Management: Publication No. EPA/530-SW-86-014. Washington, DC: US Environmental Protection Agency, 1986:1-1 ~ 5-5, R1-R3, A1-A24.
9. CDC: 2008 Case Definition; Human Immunodeficiency Virus Infection.
10. Delaney KP, Branson BM, Uniyal A, et al. Evaluation of the Performance Characteristics of 6 Rapid HIV Antibody Tests. *Clinical Infectious Diseases*. 2011; 52(2): 257-263.
11. O'Connell RJ, Merritt TM, Malia JA, et al. Performance of the OraQuick Rapid Antibody Test for Diagnosis of Human Immunodeficiency Virus Type 1 Infection in Patients with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2003; 41(5):2153-2155.
12. O'Connell RJ, Agan BK, Anderson SA, et al. Sensitivity of the Multispot HIV-1/HIV-2 Rapid Test Using Samples from Human Immunodeficiency Virus Type 1-Positive Individuals with Various Levels of Exposure to Highly Active Antiretroviral Therapy. *Journal of Clinical Microbiology*. 2006; 44(5): 1831-1833.

Advice Line

For further information, please contact your distributor, or call to one of the following Alere Product Support Care Centers:

Region	Phone	E-Mail Address
Europe & Middle East	+44 (0) 161 483 9032	EMEproductsupport@alere.com
Asia Pacific	+ (61) 7 3363 7711	APproductsupport@alere.com
Africa, Russia & CIS	+ (972) 8 9429 683	ARCISproductsupport@alere.com
Latin America	+ (57) 2 661 8797	LAprductsupport@alere.com

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Alere Ltd.

Pepper Road, Hazel Grove,
Stockport, SK7 5BW, UK

Alere Medical Co., Ltd.

357 Matsuhida, Matsudo-shi,
Chiba, 270-2214 Japan

Tel +81 47 311 5750

www.alere.com