

GOLD

Hp dry

LENCOMM S.J.

LENCOMM

Phone: +4822 8969487

Fax: +4822 8350354

e-mail: lencomm@lencomm.pl

!NEW!

For the Diagnosis of *Helicobacter pylori*



Epidemiological studies have shown a close correlation between the presence of *Helicobacter pylori* in the stomach and chronic type B gastritis, gastroduodenal ulcers and the risk of getting certain types of gastric adenocarcinomas.

GOLD HP Dry is a qualitative rapid test to detect *Helicobacter pylori* in gastric biopsies based on the activity of the enzyme urease.

The test has been designed on a polypropylene well with a urea containing cellulose support. A sample with *Helicobacter pylori* will react with the urea liberating ammonia producing a shift in pH. This shift is visualized by a change in the color of the cellulose support from yellow to purple of an indicator dye.

High Sensitivity

High Specificity

Simple

No Incubation or Reagentes Required

Fast

Easy to Read

Room Temperature Storage



www.lencomm.pl